



# [**Fr** > **Ft**] $\therefore \Delta \mathbf{F} = (\mathbf{Fr} - \mathbf{Ft}) = \underline{Positive}$ (above baseline)



# $[Fr = Ft] \therefore \Delta F = (Fr - Ft) = Zero$ (no flow detected)




 $[Fr < Ft] \therefore \Delta F = (Fr - Ft) = \underline{Negative} (below baseline)$ 

# Flash Card: Doppler Equation

# $\Delta \mathbf{F} = \mathbf{V} \operatorname{Cos} \Theta 2 \operatorname{Ft} / C$

## There is a <u>two</u> in this equation because there are two Doppler shifts.

 $\Delta F$  = Delta F = (Fr – Ft)Fr = freq receivedV = velocity of blood flow, This is what we are solving for!Ft = freq transmitted $\Theta$  = Theta =angle between the US beam and blood flow.Ft = freq transmittedC = 1.54 mm/microsecond = 1540 m/s (speed of sound in soft tissue) (memorize this)



#### PTEmasters.com

## <u>A color MAP = variance map is shown.</u>

<u>The type of flow (velocities) coded for by the</u> <u>color map are as follows:</u>

A = Black = No Doppler shift. No flow detected.

**B** = <u>Yellow</u> = Turbulent flow toward the probe

**C** = **Red** = **Laminar** flow toward the probe

- **D** = **Blue** = **Laminar** away from the probe
- **E** = **Green** = **Turbulent** away

#### <u>References:</u>

1. Understanding Ultrasound Physics by Edelman, SK, 3rd ed. Pages 314

#### **PTEmasters.com**

#### $\Delta F = V \cos \Theta 2Ft / C$

The Nyquist limit is the max Doppler shift (NOT max velocity) that can be measured. Nyquist limit =  $\frac{1}{2}$  x PRF (PRF = pulse repetition frequency)

**PRP** = pulse repetition period = inverse of PRF, thus Nyquist limit = 2 x PRP.

The Nyquist limit is a delta F = Doppler shift = a Frequency.

The Nyquist limit is also called the Nyquist frequency.

A negative Doppler shift is displayed below the baseline.


A positive Doppler shift is displayed above the baseline

Red cells moving away from the probe cause a NEGATIVE Doppler shift or Delta F.

**Red cells moving toward the probe will cause a POSITIVE Doppler shift or Delta F = (Received Freq – Transmitted Freq)** 

 $\Delta F = (\text{Received Freq} - \text{Transmitted Freq}) = (FR - FT)$ 

A <u>negative Doppler shift</u> occurs when the <u>frequency received</u> is <u>LESS</u> than the frequency transmitted. A <u>positive Doppler shift</u> occurs when the <u>frequency received</u> is <u>MORE</u> than the frequency transmitted.



References:1. Understanding Ultrasound Physics by Edelman, SK, 3rd ed. Page 2962. Mathew, JP et. al. Clinical Manual and Review of TEE, 3rd ed. McGraw-Hill 2019:13-16

# Flash Card Opportunity:

- IF Doppler angle is >0, V is underestimated.
- **Cosine 0 degrees = 1 (estimates actual velocity)**
- **Cosine 20 degrees = 0.94 (6% underestimate of V)**
- **Cosine 30 degrees = 0.87 (13% underestimation of V)**
- **Cosine 50 degrees = 0.64 (36% underestimate of V)**
- **Cosine 60 degrees = 0.50 (50% underestimate of V)**
- **Cosine 90 degrees = 0 (No Doppler shift, No Velocity detected)**

Aliasing artifact occurs if the Doppler shift has exceeded 1/2 PRF.

- <sup>1</sup>/<sub>2</sub> x PRF is known as the Nyquist limit or Nyquist frequency.
- If the PRF is increased, it means that ½ x PRF is also increased and if ½ PRF is increased it makes
- it less likely that the Doppler shift will exceed this "Nyquist limit = Nyquist frequency".
- So, if PRF is increased, <sup>1</sup>/<sub>2</sub> PRF = Nyquist limit is increased and aliasing is less likely.

<sup>1</sup>/<sub>2</sub> x PRF = Nyquist limit (Nyquist frequency) = the Doppler shift above which aliasing artifact occurs.

### 4 Ways to decrease aliasing artifact:

- **1. Decrease the depth of the sample volume**
- 2. Lower the transmitted frequency
- 3. Adjust (optimize) the baseline to focus on the flow of interest
- 4. Use continuous wave Doppler.

#### **Decreasing the transmitted frequency (Ft)** will decrease *aliasing artifact*.

I remember this in the following way (might not be perfectly correct but works for me):

IF the transmitted frequency (FT) is decreased, then the calculated Doppler shift ( $\Delta$ F) will be lower. We see this when we look

at the Doppler equation:


 $\Delta$  F = V cosine  $\Theta$  2Ft / C

<u>References:</u>

1. Understanding Ultrasound Physics by Edelman, SK, 3rd ed. Pages 308-309

If Ft is lower, then the calculated  $\Delta$  F is lower, and it is less likely to exceed the max  $\Delta$  F or Nyquist frequency = Nyquist limit =  $\frac{1}{2}$  x PRF.

## CW = Continuous Wave Doppler



[**Fr** > **Ft**]  $\therefore \Delta \mathbf{F} = (\mathbf{Fr} - \mathbf{Ft}) = \underline{Positive}$  (above baseline)

CW: TWO crystals: 1 & 2 One continuously sending (1) One continuously receiving (2) **CW Doppler involves NO damping** 

- **CW Doppler transducers have a narrow = small bandwidth.**
- **CW** Doppler transducers have a high **Q** factor = **Q**uality factor.
- **Q** factor = **Q**uality factor = **Resonant** frequency / bandwidth
- CW has NO Damping and this gives them a small band width, and increases **Q** factor.
- **CW Doppler transducers have two crystals: one continuously sending, one continuously receiving.**
- **CW Doppler can measure high velocities.**
- **CW Doppler lacks the ability to pin point the exact location of velocities**
- (No range resolution = range ambiguity).
- Aliasing does NOT occur with CW Doppler.

**References** 

1. Understanding Ultrasound Physics by Edelman, SK, 3<sup>rd</sup> ed. Page 296 . Mathew, JP et. al. Clinical Manual and Review of TEE, 3rd ed. McGraw-Hill 2019:13-16